LINFOPRESS STUDIO

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• APPLICATORS AND ACCESSORIES:

- Leg strap six sectors
- Pair of legs strap six sectors
- Arm strap six sectors
- Pair of arms strap six sectors
- Six-outlet duplication set
- Leg strap four sectors
- Pair of legs strap four sectors
- Four-outlet duplication set
- Two sectors abdominal strap
- Splitter set for two sectors abdominal strap

- CE Certified: this unit fully complies to specifications given by the European Community under 93/42/EEC e 89/336/EEC Directives.

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The technical features and sizes of the appliances described herein are subject to modification, based on technological evolution and on the constant improvement of our production.

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





SEQUENTIAL OPERATION

The treatment starts when the first sector inflates (Phase 1) until the pressure previously set by the operator is reached.

Then, the second sector inflates (Phase 2). Once Phase 2 reaches its steady state, the system reaches Phase 3 and so on until the last phase is reached.

Once all the sectors have been inflated, then they start to deflate (Phase 7).

Linfopress Studio is the last novelty in sequential pressotherapy.

This type of appliances that have been used for years for medical and physiotherapy, allows to act physiologically, wherever help for problems due to an incorrect functioning of the lymphatic or venous system is needed.

Linfopress Studio has been designed and realized to comply with professional requirements for medical and physiotherapy applications: the attractive design, compact size and ease of use makes Linfopress Studio suitable for treatments both in an ambulatory and at home.

The appliance adopts a powerful compressor which allows to obtain a 30 second inflating cycle with the possibility for the operator to adjust the treatment time from 5 to 60 min. and to adjust the pressure precisely from 0 ÷ 160 mmHg; both parameters are displayed on a large backlit screen.

The inflating cycle takes place sequentially on six different sectors of the leg and the arm strap with overlapped sector inflating system. This inflating system is particularly interesting in that it does not let any vacuum spaces build between the sectors and the massage obtained is complete, uniform and effective.

Legs and arms strap of **Linfopress Studio** are made of a special plastic material, extremely resistant and very soft, in order to avoid rubescence.

The appliance can pilot a single applicator or two applicators by the specific splitters; the leg straps allow to use a special two-sectors abdominal strap by means of an extra connection to the inguinal sector.

Linfopress Studio works in a very easy way. After applying the leg or the arm strap, start the appliance by pressing the power supply switch. The LC display will prompt the parameters the operator can set.

The treatment cycle consists in inflating the six sectors one after the other with 25 seconds' time and 5 seconds' deflating. The cycle is repeated for the time set in advance at the end of which an acoustic signal will be released. The treatment pressure in the various sectors is adjusted by the operator by acting on the knob on the front panel.

Linfopress Studio permits to act on venostasis and limphostasis, with a homogeneous and constant massage; it is useful in treating lymphoedemas and various kind of oedema; furthermore it is suitable for short run after a long physical effort in sport and to fight against the blemishes of cellulitis.

fisioline
biomedical instrumentation