

## FISIOFIELD MINI

Portable instrument for magnetotherapy with 2 outputs. It can be used at the patient's or in an ambulatory. Square wave.



12

- Magnetic field's intensity: 0-100 GAUSS for each applicator.
- Emission frequency: 1-100Hz.
- Pair of contact applicators  $\varnothing$  9 cm.

Generator size:  
220 x 156 x 100 mm.

Optional:

- Additional pair of applicators  $\varnothing$  9 cm.
- Pair of contact applicators  $\varnothing$  13 cm.
- Instrument carry-on.

