

WINNER®



Battery neuromuscular stimulator
with power supply.
4 independent stimulation channels
for band width adjustment.
Output stage: current controlled.

Wave shape:
rectangular two phase symmetrically
compensated wave.
Pulse duration:
adjustable from 50µs to 500µs.
Pulse frequency:
selectable from 1Hz to 200 Hz.
Current output:
up to 120 mA.
Graphic display.



11

Stimulation equalizer

a revolutionary device permitting to
harmoniously distribute power on all
muscle fibers for unsurpassed stimulation.

Possibility of using two different programs
on two groups of channels (function 2+2).

600 preset programs subdivided in
SPORT e FITNESS, REHAB, RELAX e TENS,
URO, BEAUTY e VASCULAR.
100 free memories to set
personalized programs.

Each program can correspond
to a sequence of as many as 5 phases.

Size: 170 x 100 x 40 mm.